

# Antibacterial agents

## Antibiotic susceptibility testing

- **Aim:** investigate the antibiotic susceptibility and resistance of pathogenic bacteria to antibiotics and determine the group of applicable antibiotics.

Group of tested antibiotics depends on

- the pathogenic bacterium and its generic resistance
- localization of infection and penetration of the drug (e.g. penetration to CSF)
- nosocomial vs. community acquired infection
- **Appropriate method should be chosen for evaluation of antibiotic susceptibility (e.g. disk diffusion method cannot be used for anaerobic bacteria).**

## Definitions

### Minimum inhibitory concentration (MIC)

the lowest concentration of an antibiotic that effectively inhibits the growth of the bacterial strain isolated from the patient (the smaller the better!)

- **MIC range**  
minimum and maximum MIC values determined in different strains of a bacterium species
- **MIC 50**  
the MIC which inhibits the growth of 50% of the strains tested
- **MIC 90**  
MIC that inhibits the growth of 90% of the strains tested

Based on the MIC range, MIC50 and MIC90, bacterial strains can be enrolled to the following categories according to their MIC values:

- **Susceptible**  
(If the MIC of the isolate is lower than a certain breakpoint value)
- **Resistant**  
(If the MIC of the isolate is higher than an other breakpoint value)
- **Intermediate**  
(If the MIC of the isolate is between the 2 breakpoint values)

Table III. Summary of breakpoint recommendations (concentration in mg/L)

	Group 1		Group 2	
	Staphylococci, streptococci, <i>M. catenulata</i> , <i>H. influenzae</i>	resistant $\geq$	Enterobacteriaceae, <i>Pseudomonas</i> spp.	resistant $\geq$
<b>5.1.1 Penicillins</b>	0.12	0.25	-	-
5.1.1.1 benzylpenicillin <sup>a</sup>	-	-	-	-
5.1.1.2 Penicillinase-resistant penicillins	-	-	-	-
nafcillin	4	8	-	-
methicillin	4	8	-	-
oxacillin	2	4	-	-
5.1.1.3 Broad-spectrum penicillins	-	-	-	-
amoxicillin	1	2	8	16
ampicillin	1	2	8	16
ox-amoxicillin	1	2	8	16
5.1.1.4 Anti-pseudomonas penicillins	-	-	-	-
piperacillin $\pm$ tazobactam	2	4	16	32
ticarcillin $\pm$ clavulanate	2	4	16	32/128 <sup>b</sup>
5.1.2 Cephalosporins, cephamycins & other $\beta$ -lactams	-	-	-	-
cefazolin	1	2	1	2
ceftriaxone	1	2	-	-
cefepime	2	4	1	2
cefuroxime	1	2	1	2
cefotaxime	1	2	1	2
cefazolin	4	8	4	8
cefuroxime	4	8	4	8
cefepime	4	8	4	8
cefoperazone	4	8	4	8
cefotaxime	1	2	1	2
ceftriaxone	1	2	1	2
cefepime	-	-	-	-
cefazolin	2	4	2/16 <sup>b</sup>	4/16 <sup>b</sup>
cefuroxime	1	2	1	2
ceftriaxone	1	2	1	2
cefuroxime iv	1	2	8	32
cefuroxime po	1	2	1	2
ceftazidime	2	4	2	4

January 2005

14

### Minimum bactericidal concentration (MBC)

the lowest concentration of an antibiotic that effectively kills at least 99.9% the bacterial isolate


#### Bactericidal agent

MBC is equal or close to the MIC

#### Bacteriostatic agent

MBC is significantly higher than the MIC

- **Serum bactericidal titre (SBT)**  
the highest dilution of the patient's serum (CSF), in which the antibiotic kills the bacterium isolated from the patient
- **Time-kill curves**  
describe the kinetics of antimicrobial effect  
show the change in the number of viable bacterium cells as a function of time  
describe the interactions between different antibiotics (antagonism, synergism, or additive interaction)



## Interactions of antibiotics

- **Synergism**  
the effect of the two drugs together is significantly higher than the sum of the effects of the two drugs used separately (penicillin + gentamicin)
- **Antagonism**  
the effect of the two drugs together is lower than the sum of the activities of the two drugs used separately (penicillin + tetracycline)
- **Indifference**  
the two drugs have additive effect when used together

## Quantitative methods Agar dilution method

- determination of minimum inhibitory concentration
- two-fold serial dilutions (decreasing concentration) of a given antibiotic in agar media
- standard inoculum of bacterium strain is inoculated onto the agar medium
- MIC is the lowest concentration of the antimicrobial drug at which no growth of bacteria can be seen
- expensive and time-consuming



Vancomycin sensitivity testing in MRSA isolates


Ann. Trop. Med. Publ. Health (2012) 5: 178-180.


## Quantitative methods Broth dilution method

- application: MIC, MBC, SBT testing
- two-fold serial dilution of the antibacterial drug in liquid medium
- standard inoculum of isolated bacterium strain is inoculated into the broth
- MIC – the lowest concentration of antibiotic that prevents the visible growth of the microorganism
- microdilution (100-200  $\mu$ l)
- macrodilution (1 ml)
- expensive and time-consuming method



## Bactericidal antibiotic





### Quantitative methods E test

- MIC
- gradient of increasing concentrations of test antibiotic is incorporated into a plastic-coated strip
- strip is placed on agar surface onto which test organism has been streaked
- antibiotic diffuses into the agar medium according to the concentration gradient
- elliptical inhibition zone
- intersection of lowest point of elliptical zone of growth inhibition and gradient strip is interpreted as minimal inhibitory concentration
- easy to use, but expensive

### Semi-quantitative methods Breakpoint method

- Antimicrobial effect is examined at only two concentrations, at the susceptible and the resistant breakpoints.
- It can be tested on solid medium and in broth containing critical concentrations of antibiotic.
- If there is no growth at either concentrations, the microbe is susceptible.
- If there is no growth at the higher concentration, but there is growth at the lower concentration, the microbe is intermediate.
- If growth can be seen at both concentrations, the microbe is resistant.

Before inoculation	<table border="1" style="display: inline-table;"> <tr><td>8</td><td>16</td></tr> <tr><td>8</td><td>16</td></tr> <tr><td>4</td><td>8</td></tr> </table>	8	16	8	16	4	8	Ampicillin Cephalothin Gentamicin	
8	16								
8	16								
4	8								
After inoculation	<table border="1" style="display: inline-table;"> <tr><td>8</td><td>16</td></tr> <tr><td>8</td><td>16</td></tr> <tr><td>4</td><td>8</td></tr> </table>	8	16	8	16	4	8	Ampicillin <b>R</b> Cephalothin <b>MS</b> Gentamicin <b>S</b>	
8	16								
8	16								
4	8								

### Semi-quantitative methods Resistance screening

- Similar to the breakpoint method but only one agar plate is used with the higher breakpoint concentration.
- If bacterial growth can be seen, the microbe is suspected to be resistant, but this should be confirmed by MIC tests.
- Is useful in heteroresistance detection and to identify aminoglycoside resistance of enterococci.

### Semi-quantitative methods Disk diffusion test (Kirby-Bauer test)


- Antibiotic-impregnated disks are placed onto the inoculated agar surface
- Incubation at 35-37°C for 18-24h.
- Antibiotics diffuse into the medium.
- If the bacterium is susceptible, zone of inhibition will be formed around the disks.
- The size of zones depends on:
  - Concentration of the antibiotic.
  - The effect of the antibiotic against microbe.
  - The diffusion ability of the antibiotic in the medium.
- The size of the zone of inhibition is compared to standards to determine the susceptibility of the microbe to the drug.
- Not useful in anaerobic infections, slowly growing microbes and fungi, but cheap and easy to perform.

### Detection of β-lactamases

- It is important if it is more reliable than culturing methods (*N. gonorrhoeae*).
- The absence of detectable β-lactamase activity does not mean that the strain is susceptible to β-lactams.
- Different methods:
  - Nitrocefin test: N. is a cephalosporin and if the lactam ring is hydrolyzed by the β-lactamase of the bacterium, colour change from yellow to red can be seen. *Haemophilus*, *Neisseria*, *Moraxella*.
  - Iodometric test: if the β-lactam ring is hydrolyzed, reduction of iodine can be detected with starch.
  - Acidimetric test: in case of β-lactam hydrolysis, pH decreases because of carboxyl groups, which can be detected with an indicator

## Detection of extended spectrum B-lactamases

- ESBL: is able to hydrolyze 3rd generation cephalosporins.
- ESBLs can be inhibited by  $\beta$ -lactamase inhibitors in vitro. So the cephalosporins will be effective in the presence of the inhibitors in vitro (but not in vivo).




### Disk diffusion method

- uses a disk only with cephalosporin and an other with cephalosporin and clavulanic acid on the same inoculated solid medium
- test is positive, if the difference between the inhibitory zones is more than 5 mm

### Double disk method

- uses a disk impregnated with cephalosporin, and an other with clavulanic acid and the disks are 2 cms apart
- the inhibitory zone stretches toward the disk containing clavulanic acid.

- In case of **E-test**, one end is impregnated only with cephalosporin, and the other with cephalosporin and clavulanic acid.



## PBP2a detection of methicillin resistant Staphylococcus aureus

- Using monoclonal antibodies in latex agglutination test.



Direct detection of resistance genes can be done by PCR